



# ALGALTOXKIT F

## Test procedure



1

### PREPARATION OF ALGAL CULTURING MEDIUM

- VOLUMETRIC FLASK (1 liter)
- VIALS WITH NUTRIENT STOCK SOLUTIONS A (2 vials), B, C, D
- DISTILLED (or deionized) WATER



2

TRANSFER 10 ML FROM ONE OF  
THE TWO "NUTRIENT STOCK A" VIALS  
IN ± 800 ML DISTILLED WATER  
IN THE 1 LITER VOLUMETRIC FLASK



3

TRANSFER 1 ML FROM THE  
NUTRIENT STOCK VIALS  
B, C AND D INTO THE 1 LITER  
VOLUMETRIC FLASK.



4

- FILL THE FLASK TO THE 1 LITER MARK WITH DEIONIZED WATER
- STOPPER THE FLASK AND SHAKE THOROUGHLY TO HOMOGENIZE THE CONTENTS
- AERATE THE ALGAL CULTURING MEDIUM FOR AT LEAST 30 MINUTES



5

ADJUST THE pH  
(if necessary)  
TO  $8,1 \pm 0,2$   
(with either 1 M HCl  
or 1 M NaOH)



6


### DE-IMMOBILIZATION OF THE ALGAE

TAKE ONE TUBE CONTAINING ALGAL BEADS AND POUR OUT THE LIQUID  
TAKE CARE NOT TO ELIMINATE ANY OF THE ALGAL BEADS DURING THE PROCESS !!



7

OPEN THE VIAL "MATRIX DISSOLVING MEDIUM" AND TRANSFER 5 ML  
TO THE TUBE WITH ALGAL BEADS



8

CAP THE TUBE AND SHAKE VIGOROUSLY TO DISSOLVE THE ALGINATE MATRIX OF THE ALGAL BEADS, PREFERABLY WITH THE AID OF A VORTEX SHAKER



9

CONTINUE THE SHAKING UNTIL THE ALGAL BEADS ARE TOTALLY DISSOLVED



10

CENTRIFUGE THE TUBE FOR 10 MINUTES AT 3000 RPM IN A CONVENTIONAL LAB CENTRIFUGE




11

POUR OUT THE SUPERNATANT FROM THE TUBE



12

- ADD 10 ML DISTILLED WATER TO THE TUBE
- CAP AND SHAKE THE TUBE TO RESUSPEND THE ALGAE



13


- CENTRIFUGE THE TUBE AGAIN AT 3000 RPM FOR 10 MINUTES  
AND THEN POUR OUT THE SUPERNATANT



14

- ADD 10 ML ALGAL CULTURING MEDIUM TO THE TUBE

- CAP THE TUBE AND SHAKE TO RESUSPEND THE ALGAE



15

### PREPARATION OF CONCENTRATED ALGAL INOCULUM


TRANSFER THE ALGAL SUSPENSION FROM THE TUBE INTO A 25 ML CALIBRATED FLASK






## 16

- ADD ALGAL CULTURING MEDIUM TO THE 25 ML MARK OF THE FLASK
- STOPPER THE FLASK AND SHAKE TO HOMOGENIZE THE ALGAL SUSPENSION




## 17

- PUT 25 ML ALGAL CULTURING MEDIUM IN THE CALIBRATION LONG CELL AND CLOSE THE CELL WITH THE LID
- PLACE THE CELL IN THE SPECTROPHOTOMETER




18

ZERO-CALIBRATE THE INSTRUMENT AT A WAVELENGTH OF 670 NM



19

TRANSFER THE 25 ML ALGAL SUSPENSION INTO THE ALGAL STOCK CELL




20

CLOSE THE ALGAL STOCK CELL WITH THE LID  
AND SHAKE TO DISTRIBUTE THE ALGAE EVENLY



21


PUT THE ALGAL STOCK CELL IN THE SPECTROPHOTOMETER  
AND READ THE OPTICAL DENSITY (**OD1**) AFTER 10 SECONDS



22


- TAKE THE **OD/N**  
 (optical density/algal number)  
 SHEET

- WITH THE AID OF THE REGRESSION FORMULA CALCULATE THE NUMBER OF ALGAE **N1** CORRESPONDING TO THE MEASURED **OD1** IN THE ALGAL STOCK CELL
- WITH **N2** =  $1 \cdot 10^6$  ALGAE/ML, CALCULATE FROM THE **N1/N2** RATIO THE DILUTION FACTOR NEEDED TO REACH **OD2** (corresponding to  $1 \cdot 10^6$  algae/ml)



23

POUR THE 25 ML ALGAL  
 SUSPENSION FROM THE  
 ALGAL STOCK CELL INTO  
 A 100 ML FLASK



24


ADD THE CALCULATED VOLUME OF ALGAL CULTURING MEDIUM  
TO THE FLASK, TO MAKE UP A SUSPENSION OF  $1.10^6$  ALGAL CELLS / ML



25


**PREPARATION OF THE  
TOXICANT DILUTION SERIES**

TAKE SIX 200 ML CALIBRATED  
FLASKS AND LABEL THEM FROM  
C0 TO C5



26

TO ELIMINATE TURBIDITY,  
SAMPLES MUST BE  
FILTERED BEFORE TESTING  
(e.g. over a membrane filter  
of 0.45  $\mu\text{m}$ ),




27

PUT 100 ML ALGAL CULTURING MEDIUM  
IN THE 200 ML FLASKS C0, C2, C3, C4 AND C5




28

FILL FLASK C1  
TO THE 200 ML MARK  
WITH THE FILTERED SAMPLE




29

- ADD 2 ML "NUTRIENT STOCK SOLUTION A" AND 0.2 ML OF NUTRIENT STOCK SOLUTIONS B, C AND D TO FLASK C1
- STOPPER THE FLASK AND SHAKE TO MIX THE CONTENTS



## 30


- POUR 100 ML SAMPLE FROM FLASK C1 INTO A GRADUATED CYLINDER AND TRANSFER THIS VOLUME INTO FLASK C2 TO MAKE THE FIRST 1:1 DILUTION
- STOPPER FLASK C2 AND SHAKE TO HOMOGENIZE THE CONTENTS



## 31

REPEAT THE FORMER DILUTION PROCEDURE FOR THE OTHER FLASKS (i.e.. 100 ml from C2 to C3, etc.)






32

REMOVE AND DISCARD  
100 ML SOLUTION FROM  
FLASK C5 TO HAVE 100 ML  
SOLUTIONS IN EACH FLASK




33

- TAKE THE FLASK CONTAINING THE  $1.10^6$ /ML ALGAL SUSPENSION AND SHAKE IT GENTLY
- ADD 1 ML ALGAL SUSPENSION TO EACH OF THE 6 FLASKS C0 TO C5, IN ORDER TO OBTAIN AN INITIAL ALGAL CONCENTRATION OF  $1.10^4$  CELLS/ML IN EACH FLASK



34

TRANSFER OF THE  
ALGAE-TOXICANT  
DILUTIONS INTO  
THE TEST VIALS



35

- LABEL ALL THE LONG CELLS ON THEIR LID (3 replicates per dilution)
- AFTER THOROUGH SHAKING, TRANSFER 25 ML ALGAE-TOXICANT DILUTION FROM EACH FLASK INTO A GRADUATED CYLINDER, FOR FURTHER TRANSFER INTO THE CORRESPONDING LONG CELL (3 replicates per toxicant dilution)



## 36

- REDISTRIBUTE THE LONG CELLS IN THE HOLDING TRAY IN A RANDOM WAY
- LIFT UP THE LIDS OF THE CELLS A LITTLE AT ONE SIDE, AND SLIDE THE PLASTIC STRIP OVER THE OPEN PART OF THE LONG CELLS TO KEEP THEM SLIGHTLY OPEN DURING THE INCUBATION PERIOD



## 37

INCUBATE THE HOLDING TRAY WITH THE LONG CELLS FOR 72h IN AN INCUBATOR AT  $23^{\circ}\text{C} \pm 2^{\circ}\text{C}$ , WITH CONTINUOUS ILLUMINATION:

- SIDEWAY ILLUMINATION = 10000 LUX
- OR BOTTOM ILLUMINATION = 3000-4000 LUX



38


### SCORING OF THE RESULTS

THE **OD** OF THE ALGAL SUSPENSIONS SHALL BE MEASURED EACH DAY DURING THE 3 DAYS OF THE TEST, I.E. AFTER 24h, 48h AND 72h EXPOSURE TO THE TOXICANT



39

ZERO-CALIBRATE THE SPECTROPHOTOMETER PRIOR TO THE DAILY MEASUREMENT OF THE **OD** IN THE LONG CELLS,



40

IMMEDIATELY BEFORE MEASURING THE **OD** IN A LONG CELL, CLOSE THE CELL, TURN IT UPSIDE DOWN AND SHAKE GENTLY TO RESUSPEND THE ALGAE EVENLY



41

- SCORE THE DAILY **OD** RESULT OF EACH LONG CELL ON THE "RESULTS SHEET"
- PERFORM THE DATA TREATMENT OF THE RESULTS WITH AN APPROPRIATE PROGRAMME