

■ ■ ■ ■ ■ Lumitester PD-30 + LuciPac Pen / LuciPac Pen AQUA LuciPac Pen-LS at a glance

Lumitester PD-30

- NEW** ■ Temperature compensation
- NEW** ■ Self-diagnosis
- NEW** ■ 8 languages
- NEW** ■ Data-analysis software
 - Reliable monitoring of liquids, and dry and moist surfaces
 - Patented measurement of ATP and AMP increases the sensitivity
 - 2,000 memory spaces
 - 400 modes
 - Measurement in RLUs – relative light units
 - Smallest device in the world
 - Lightest device in the world
 - Ready for use in 8 seconds
 - Auto-calibration
 - For auditing, the device can be sent to us for recalibration and certification
 - Data transmission via USB to PC
 - Documentation of the results
 - Optional graphic analysis of the data
 - Excel-based software is easy to operate
 - Software in English facilitates international use
 - Definition of your own limit values to monitor individual checking samples
 - 2 year guarantee
 - CE mark and produced to ISO 9001
 - Meets DIN 10214: 2009-12
 - Supplied with self-supporting protective cover

LuciPac Pen / LuciPac Pen AQUA / LuciPac Pen LS (Long Swab)

- **LuciPac Pen:** For dry and moist surfaces
 - **LuciPac Pen AQUA:** For liquids and water
 - **LuciPac Pen LS:** For endoscopes and hard-to-access locations
 - Heat-sealed and unbreakable swab
 - Swabs can be removed individually from resealable packaging
 - Detergent-tolerant enzymes – low inhibition from detergent residues
 - High stability of swab at room temperature
 - Greatest sensitivity due to patented simultaneous measurement of AMP and ATP
 - Shelf life: 15 months at 2 – 8°C; 14 days at 25°C; 5 days at 30°C.
- The product has a shelf life of up to 30 days at temperatures of 20°C or less.
- If the swabs are cooled again, the labelled shelf life of 15 months then reapplies!



LuciPac Pen-LS (Long Swab)

